Bluephase[®] Meter II

EN Instructions for Use

DE Gebrauchsinformation

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SV Bruksanvisning

DA Brugsanvisning

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NO Bruksanvisning

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For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Austria

Ivoclar Vivadent AG

FL-9494 Schaan / Liechtenstein

In USA distributed by

Ivoclar Vivadent Inc., 175 Pineview Drive,
Amherst, NY 14228, USA

CE

ivoclar[®]
vivadent[®]
clinical

Dear Customer

Thank you for having purchased the Bluephase Meter II from Ivoclar Vivadent. The light intensity of the polymerization light used is the decisive factor when it comes to the sufficient curing of restorations. In order to achieve adequate polymerization and thus durable high-quality composite restorations, it is recommended to check the light intensity of the polymerization lights regularly with a radiometer.

These Instructions for Use will help you to correctly start-up the apparatus, to make full use of its capabilities and to ensure a long service life. Should you have any further questions, please do not hesitate to contact us (see addresses on last page).

Intended use

The Bluephase Meter II is a radiometer for measuring energy-rich blue light in the wavelength range of 380–550 nm in polymerization lights intended for dental use. The Bluephase Meter II is suitable for measuring the light intensity (irradiance) in [mW/cm²] as well as the light power (radiant flux) in [mW]. The intended use includes the observation of the notes and regulations in these Instructions for Use. Position of use – lying on the table.

Signs and symbols

Contraindication

Symbols on the device

Observe Instructions for Use

The apparatus must not be disposed in the normal domestic waste. Please visit your local Ivoclar Vivadent website for detailed information about the disposal of the unit.

Safety notes

The Bluephase Meter II is a laboratory apparatus which is subject to IEC 61010-1 (EN 61010-1). It complies with the relevant EU regulations. The apparatus has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use have to be observed.

NRTL (Nationally Recognized Testing Laboratory)

Warning

- Protect the measuring device from heavy shock or falling down.
- Cleaning agents and other liquids must not enter the apparatus, since they may damage it.
- Do not expose the apparatus to sunlight.
- Both the storage and use of other electronic devices in the immediate vicinity of the apparatus may impair proper operation.

- Portable and mobile high-frequency communication devices may affect the apparatus. The use of mobile phones during operation is not allowed.
- Repair work may only be carried out by a certified Ivoclar Vivadent Service Center.
- For dental use only! The Bluephase Meter II is suitable for measuring energy-rich blue light in the wavelength range of 380–550 nm and not white light.

Start-up

To start the apparatus, open the battery case at the rear of the housing and insert the provided batteries (3 x type LR6/AA/1.5V). Do not use any rechargeable batteries. The unit is now ready for operation.

Operation

For the correct measuring with the Bluephase Meter II, the design of the polymerization light must be considered. For polymerization lights with the light source positioned directly at the light emission window of the handpiece (e.g. LED at the front of the light head), only the light power (radiant flux) in [mW] can be measured.

Measuring the light intensity (irradiance) in [mW/cm²]

- 

1.) Define the diameter of the circular light probe using the integrated template on the rear of the device.
- 

2.) Set the defined diameter between << Ø 6 – 12 mm >> using the (-) or (+) buttons.
- 

3.) Position the light probe directly on the centring gauge and flush on the sensor surface.

- 4.) Switch on the polymerization light. The Bluephase Meter II is activated automatically. The previously set diameter is displayed, providing visual feedback.
- 5.) The measurement result of the light intensity is displayed in milliwatt per square centimetre [mW/cm^2]. **Note:** If the light intensity is below $300 \text{ mW}/\text{cm}^2$, "MIN" is shown in the display.
- 6.) Turn off the polymerization light. The Bluephase Meter II turns off automatically.

Measuring the light power (radiant flux) in [mW]

- 1.)  Press the (+) button until the milliwatt program is selected and << ----mW >> is shown in the display.
- 2.)  Corresponds with Steps 3.) to 6.), the measurement result of the light performance is displayed in milliwatt [mW].

To convert [mW] to [mW/cm^2], the measured light intensity must be divided by the active light emission surface [cm^2].

Example: 670 mW divided by 0.608 cm^2 results in $1,102 \text{ mW}/\text{cm}^2$.

Acoustic signal

If any button is pressed, an acoustic signal sounds. By pressing the (+) and (-) buttons simultaneously and for a longer period, the signal tone can be activated or deactivated.

What if ...

- ... the measurement results correspond with the data provided by the manufacturer?
The polymerization light achieves the stipulated light intensity and can therefore be used according to the polymerization recommendations of the respective manufacturer or stipulations of the corresponding composite manufacturer.
- ... the Bluephase Meter II shows a value below $500 \text{ mW}/\text{cm}^2$?
 $500 \text{ mW}/\text{cm}^2$ is the absolute minimum of light intensity for polymerization lights. If the light intensity is lower, it cannot be ensured that the composite material is completely polymerized. Therefore, it is recommended to replace a curing light with a light intensity of less than $500 \text{ mW}/\text{cm}^2$.

- ... **the measurement results deviate strongly from the data provided by the manufacturer?**

The following procedure is recommended:

- Check the selected curing program, i.e. check if the correct diameter of the light probe is set or the milliwatt program is selected.
- Clean the sensor of the Bluephase Meter II (see Maintenance and Cleaning).
- Clean the light probe with a disinfectant and a soft cloth. Carefully remove any possible material residue, such as composite material, from the light probe with a fingernail or a plastic spatula, for instance. Do not use sharp or pointed instruments, as they may scratch the surface of the light probe and thus reduce light transmission.
- Check the light probe for damage. If certain segments appear black once the light probe is removed from the handpiece and held against a light source, then the glass fibres are broken. Due to the reduced light transmission, the light probe has to be replaced with a new one.

If these measures do not improve the results, the curing times have to be prolonged in order to achieve sufficient polymerization. The following guideline applies:

If the light intensity is reduced by half, the curing time has to be doubled (e.g. the measured light intensity is 500 mW/cm² while the data stipulated by the manufacturer is 1,000 mW/cm² with 10 seconds of curing time; in this case, the curing time has to be doubled to 20 seconds).

- ...  **is shown on the display?**

A low residual capacity of the battery is signalled with a blinking battery symbol on the display. The batteries have to be replaced as soon as possible.

- ... **the display does not show anything?**

- Check or replace the batteries.
- If the display does not show anything even when exposed to the light of a polymerization device, send the Bluephase Meter II to a certified Ivoclar Vivadent Service Center for repair.

Change of batteries

Batteries are consumables and must be replaced regularly.

For this purpose, remove the cover of the battery case at the rear of the unit and replace the batteries. Subsequently,

 remount the cover of the battery case. Dispose of unserviceable batteries according to the national rules and regulations.

Maintenance and cleaning

Clean the Bluephase Meter II with a soft cloth and a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soap solution.

- Do not spray disinfecting solution directly on the apparatus and do not hold it under running water in order to prevent liquids from entering the device.
- Do not sterilize the Bluephase Meter II in the autoclave.

Warranty

The warranty for the Bluephase Meter II is three years from the date of purchase. Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases.

This is particularly valid for:

- damage resulting from improper handling;
- damage to components resulting from wear under standard operating conditions (e.g. sensor);
- damage resulting from external influences, e.g. blows, drop to the floor;
- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centers.

In case of a claim under warranty, the apparatus must be sent to the dealer or directly to Ivoclar Vivadent in its original packaging together with the purchase document.

Procedure in case of repair

Repair work may only be carried out by a certified Ivoclar Vivadent Service Center. In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on last page). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

Product specifications

Delivery form

- 667124 Bluephase Meter II
- 667237 Instructions for Use Bluephase Meter II (WE3)
- 3 x batteries (type LR6/AA/1.5V)

Technical data

Operating voltage	4.5 VDC
Battery operation	3 x type LR6/AA/1.5V

Measuring range

Wavelength range	380–550 nm
------------------	------------

Light intensity	300 – 12,000 mW/cm ²
-----------------	---------------------------------

Polymerization lights with circular light probes	Ø 6–12 mm, measured value in [mW/cm ²] and [mW]
--------------------------------------------------	-------------------------------------------------------------

Polymerization units with LED attached to the tip	Ø 5–13 mm, measured value in [mW]
---------------------------------------------------	-----------------------------------

Measuring accuracy	± 10%
--------------------	-------

Operating conditions

Temperatur	+10 °C to +35 °C (+50 °F to +95 °F)
------------	-------------------------------------

Relative humidity	30% to 75%
-------------------	------------

Ambient pressure	700 hPa to 1,060 hPa
------------------	----------------------

Dimensions :	L = 126 mm; W = 80 mm; H = 41 mm
--------------	----------------------------------

Weight (incl. batteries)	170 g
--------------------------	-------

Transportation and storage conditions

- Temperature	-20 °C to +60 °C (-4 °F to +140 °F)
---------------	-------------------------------------

- Relative humidity	10% to 75%
---------------------	------------

- Ambient pressure	500 hPa to 1,060 hPa
--------------------	----------------------

- Store the apparatus in closed, roofed rooms.

- Protect the apparatus from severe jarring.

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Ivoclar Vivadent Limited

Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE
United Kingdom
Tel. +44 116 284 7880
Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us

The logo features a series of seven dots of varying sizes arranged in a slight arc above the text. The text 'ivoclar' is in a lowercase, sans-serif font, followed by 'vivadent' in a larger, bold, lowercase sans-serif font. A registered trademark symbol (®) is located to the upper right of 'vivadent'. Below this, the word 'clinical' is written in a smaller, lowercase, italicized sans-serif font.

ivoclar
vivadent®
clinical